

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Asiga Max	385 nm	6.1 mW/cm ²	Composer 1.2.11	Detax_Medicalprint mould COLOR_5
	Asiga Max 2	385 nm	7.0 mW/cm ²	Composer 2.0.8	Detax_Medicalprint mould COLOR_5
	Asiga Pico 2	385 nm	20 mW/cm ²	Composer 1.2.11	Detax_Medicalprint mould COLOR_5
	Asiga PRO 2	385 nm	5.7 mW/cm ²	Composer 1.2.11	Detax_Medicalprint mould COLOR_5
	Asiga PRO 4K	385 nm	7.0 mW/cm ²	Composer 1.2.11	Detax_Medicalprint mould COLOR_5
	Asiga Ultra	385 nm	6.6 mW/cm ²	Composer 2.0.8	Detax_Medicalprint mould COLOR_5
	Microlay Versus	385 nm	4.3 mW/cm ²	Microform 1.0.3.7	DETAX Medicalprint Mould COLOR 100 microns v5.3

*The set of parameters includes all relevant material- and printer specific information

The placeholder **COLOR** shall be replaced by **clear**, **clear 2.0**, **rose**, **rose 2.0** or **rose-orange 2.0**

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Miicraft Alpha	385 nm	6.0 mW/cm ²	Utility 6.4.4	DETAX Medicalprint mould clear 2.0
	Miicraft Prime Series	385 nm	6.2 mW/cm ²	Utility 6.3.0	Detax_Medicalprint mould COLOR_100_PH
	Miicraft Ultra Series	385 nm	5.7 mW/cm ²	Utility 6.3.0	Detax_Medicalprint mould COLOR_100
	Rapidshape HA20+/ HA30+/HA40+	385 nm	2.0 mW/cm ²	Netfabb 2020	MedicalPrint mould COLOR
	Rapidshape HA70+/HA90+	385 nm	5.0 mW/cm ²	Netfabb 2020	MedicalPrint mould COLOR
	Way2Produc- tion SolFlex Series	385 nm	8.0 mW/cm ²	Netfabb 2020	MedicalPrint mould COLOR

*The set of parameters includes all relevant material- and printer specific information
The placeholder **COLOR** shall be replaced by **clear**, **clear 2.0**, **rose**, **rose 2.0** or **rose-orange 2.0**

Cleaning Equipment

Cleaning unit	
Manufacturer, Model	Cleaning process**
Rapidshape RS wash	Use the following settings: MedicalPrint mould COLOR Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

The placeholder **COLOR shall be replaced by **clear**, **clear 2.0**, **rose**, **rose 2.0** or **rose-orange 2.0**



detax GmbH

Carl-Zeiss-Str. 4 • 76275 Ettlingen • Germany
+49 7243 510 0 • post@detax.com • detax.com

07/2025

Light curing Equipment

Light Curing unit

Manufacturer, Model Curing process***

Dentalfarm Photopol	2 × 6 min, 120 % + N ₂ , turn around components after 6 min
------------------------	--

NK Optik Otoflash G171	2 × 2000 flashes under inert gas, turn around components after 2000 flashes
---------------------------	---

NK Optik Otoflash 250/500	4000 flashes under inert gas @15 Hz
------------------------------	-------------------------------------

Rapidshape RS cure	Use the following settings: MedicalPrint mould COLOR
--------------------	---

Rapidshape RS cure XL	Use the following settings: MedicalPrint mould COLOR
--------------------------	---

***The placeholder **COLOR** shall be replaced by **clear**, **clear 2.0**, **rose**, **rose 2.0** or **rose-orange 2.0**



detax GmbH

Carl-Zeiss-Str. 4 • 76275 Ettlingen • Germany
+49 7243 510 0 • post@detax.com • detax.com

07/2025

Mechanical finishing

Polishing unit

Manufacturer, Model	Polishing process
Otec ECO-Maxi wet	Step 1: Media: DZS 10/10 Ceramic grinding wheel, Compound: SC 15 Compound, Runtime: 90 min, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 4 l/h
	Step 2: Media: KM 10 plastic abrasive media (70 %), PM 10 plastic abrasive media (30 %), Compound: SC 15 Compound, Runtime: 90 min, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 4 l/h
	Step 3: Media: ZSP 4/5 Porcelain polishing body, Runtime: 60 min, Splitting system: Wet splitting, Speed: 260 U/min, Water flow: 10 l/h
Otec CF 3x18	Step 1: Media: DZS 10/10 Ceramic grinding wheel, Compound: SC 15 Compound, Runtime: 90 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 8 l/h
	Step 2: Media: KM 10 plastic abrasive media (70 %), PM 10 plastic abrasive media (30 %), Compound: SC 15 Compound, Runtime: 90 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 8 l/h
	Step 3: Media: ZSP 4/5 Porcelain polishing body, Runtime: 60 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 240 U/min, Water flow: 15 l/h



detax GmbH

Carl-Zeiss-Str. 4 • 76275 Ettlingen • Germany
+49 7243 510 0 • post@detax.com • detax.com

07/2025