detax

3D Medicalprint® Material

nobreak

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Intended use

Resin material for the production of earmolds and ITE shells in 3D printing processes

Indication

Earmolds, in-ear monitoring and hearing aid shells for connection to a hearing system

Patient Target Group

Persons for whom an earmold is to be created.

Intended Users

Earmold laboratory, ENT specialist, hearing aid acoustician

Suitable for the following dlp printers/cleaning/post-exposure

see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Lacquer surface or polish surface mechanically.

Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- The processing instructions and precautionary measures must be observed. Otherwise, irreparable damage to the hearing organ or eardrum may result.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly closed, carefully close immediately after each use.
- Read and understand the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- nobreak is to be stored dry (at 15 °C – 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of nobreak may cause allergic reactions in predisposed persons. In such cases refrain from using the product. nobreak only insert in completely polymerised state.

Adverse effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage

15 °C
59 °F

28 °C
82 °F

Application

At 23 °C ± 2 °C

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

① Construction process

Generation of a Print Job complying with machine and material parameters

② Post-processing

If possible, post-processing should commence immediately following the construction process. After raising the platform, a drip time of approx. 10 minutes is recommended.

③ Cleaning

see „Annex 1, Cleaning equipment“

④ Post-exposure

see „Annex 1, Light curing equipment“

⑤ Surface processing

Varnishing (e.g. with Medicalprint® coat / soft coat) or mechanical finishing* (polishing) of the surface

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detax GmbH