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Light-curing resin based on (meth)acrylate, for cast parts in dental precision casting, burn-out for DLP printers with UV-LED 385 nm

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

Processing

- ▶ Phosphate-bonded investment material can be used for fabrication of the casting mould. A preheating temperature of approx. 800 °C, is to be maintained for approx. 30–45 minutes, depending on the size of the casting mould, regardless of the final casting temperature.
- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Processing temperature 23 °C ± 2 °C.
- ▶ A correction or repair can be done with easyform gel LC after printing.

Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® CAST 2.0** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

Notes

- ▶ DETAX does not accept liability for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

Storage

- ▶ **FREEPRINT® CAST 2.0** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Area of application:

Fabrication of residue-free combustible cast objects for precision casting technology employing the DLP process.

Processing:

at 23 °C ± 2 °C

Storage

Ordering information:

FREEPRINT® CAST 2.0 385
500 g, red **02548**
1.000 g, red **02632**

FREEPRINT® MODEL 2.0 385

1.000 g
sand **02128**
light-grey **02099**
grey **02177**
caramel **02850**
white **02148**

500 g
caramel **02844**

FREEPRINT® GINGIVA 385

500 g **02820**
1.000 g **02843**

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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DETAX GmbH