FREEPRINT®

3D PRINTING MATERIAL

DENTAL MED RESIN

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**REMOVABLE
DENTURE BASES**

DETAX

de	Gebrauchsanweisung	3
en	Instructions for use	6
fr	Mode d'emploi	9
es	Modo de empleo	12
it	Istruzioni per l'uso	15
da	Brugsanvisning	18
el	Οδηγίες χρήσης	21
hr	Upute za uporabu	24
lv	Lietošanas instrukcija	27
lt	Naudojimo instrukcija	30
nk	Gebruiksaanwijzing	33
no	Brugsanvisning	36
pl	Instrukcja stosowania	39

pt	Manual de instruções	42
ro	Instrucțiuni de utilizare	45
sv	Bruksanvisning	48
sk	Návod na použitie	51
sl	Navodila za uporabo	54
cs	Návod k použití	57
hu	Használati útmutató	60
ru	Инструкция по применению	63
tr	Kullanım kılavuzu	66
fi	Käyttöohje	69
et	Kasutusjuhend	72
bg	Инструкция за употреба	75

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Denture bases

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist, dental technician

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

PROCESSING

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 8).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Polish surface mechanically.
- ▶ Processing temperature $23\text{ °C} \pm 2\text{ °C}$.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.
- ▶ MMA free.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

CLEANING

By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

By the patient

For daily cleaning of the prosthesis, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards. Denture cleaning tablets can be used following the manufacturer's instructions.

Note

- ▶ Do not use toothpaste; the abrasive particles are coarse and scratch the denture.
- ▶ Do not use any acids or solvents.

NOTES

- ▶ DETAX does not accept liability for any damage caused by misuse of the impression material.
- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.
- ▶ Read and understand the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

STORAGE

- ▶ **FREEPRINT® DENTURE IMPACT** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

CONTRAINDICATION

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **FREEPRINT® DENTURE IMPACT** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® DENTURE IMPACT** only insert intraorally in completely polymerised state.

ADVERSE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage:

Processing:

at 23 °C ± 2 °C

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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DETAX GmbH