

# Annex 1







## Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer	Light source	Light intensity	Operation Software	Parameter data set*
	Asiga Max	385 nm	6.1 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint ortho 385_5
	Asiga Max 2	385 nm	7.0 mW/cm <sup>2</sup>	Composer 2.0.8	Detax_Freeprint ortho 385_5
	Asiga Pico 2	385 nm	20 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint ortho 385_5
	Asiga PRO 2	385 nm	5.7 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint ortho 385_5
	Asia PRO 4K	385 nm	7.0 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint ortho 385_5
	Asiga Ultra	385 nm	6.6 mW/cm <sup>2</sup>	Composer 2.0.8	Detax_Freeprint ortho 385_5

\*The set of parameters includes all relevant material- and printer specific information

# Annex 1

## Digital light Processing (DLP) printer, operation software and parameter






Model (Picture)	Printer Manufacturer	Light source	Light intensity	Operation Software	Parameter data set*
	Ivoclar PrograPrint PR5	388 nm	16 mW/cm <sup>2</sup>	PrograPrint CAM 1.1.10.1	Detax FREEPRINT ortho
	Microalay Versus	385 nm	4.3 mW/cm <sup>2</sup>	Microform 1.0.3.7	DETAX Freeprint Ortho 385 100 microns v5.3
	Miicraft Alpha	385 nm	6.0 mW/cm <sup>2</sup>	Utility 6.4.4	DETAX freeprint ortho
	Miicraft Prime Series	385 nm	6.2 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint ortho 385_100_PH
	Miicraft Ultra Series	385 nm	5.7 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint ortho UV_100
	Phrozen Sonic 4K 2022**	405 nm	2.0 mW/cm <sup>2</sup>	Phrozen DS Slicer	DETAX Freeprint ortho 405

\*The set of parameters includes all relevant material- and printer specific information

\*\* The printers only works with the Freeprint ortho 405.

# Annex 1

## Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer	Light source	Light intensity	Operation Software	Parameter data set*
	Rapidshape D10+ / D20+/ D30+ / D40+	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint- ortho uv
	Rapidshape D70+ / D90+	385 nm	5.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint- ortho uv
	Rapidshape Pro20 (Dental)/ Pro30 (Dental)	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint- ortho uv
	Shining 3D AccuFab L4D**	405 nm	2.0 mW/cm <sup>2</sup>	AccuWare	DETAX FREEPRINT ortho 405
	Way2Production SolFlex Series	385 nm	8.0 mW/cm <sup>2</sup>	Netfabb 2020	Freeprint ortho 385

\*The set of parameters includes all relevant material- and printer specific information

\*\* The printers only works with the Freeprint ortho 405.

# Cleaning Equipment

---

## Cleaning unit

Manufacturer, Model	Cleaning process
---------------------	------------------

---

<b>Ivoclar PrograPrint Clean</b>	Clean the parts with isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.  The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes.  Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
--------------------------------------	---

---

<b>Rapidshape RS Wash</b>	Use the following settings: DETAX Freeprint-ortho uv.  Prior to post-exposure, check the openings, cavities and gap areas for residues.  Then blow off with compressed air.
-------------------------------	---

---

<b>Ultrasonic bath Bandelin Sonorex</b>	Clean the parts with isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.  The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
---	--

---



detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany  
+49 7243 510 0 · post@detax.com · detax.com

04/2026

## Light curing Equipment

---

### Light Curing unit

#### Manufacturer, Model

#### Curing process

---

#### Asiga Cure

Use the following settings: Detax\_Freeprint ortho 385\_6

---

#### Ivoclar PrograPrint Cure

Post curing A: Wavelength = 405 nm; Intensity = 100 %;  
Duration = 120 s  
Post curing B: Wavelength = 460 nm; Intensity = 100 %;  
Duration = 120 s

---

#### NK Optik Otoflash G171

2 × 2000 flashes under inert gas, turn around components after 2000 flashes

---

#### NK Optik Otoflash 250/500

4000 flashes under inert gas @15 Hz

---

#### Rapidshape RS Cure

Use the following settings: DETAX Freeprint-ortho uv

---

#### Rapidshape RS cure XL

Use the following settings: DETAX Freeprint-ortho uv

---



detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany  
+49 7243 510 0 · post@detax.com · detax.com

04/2026