detax

3D Freeprint® Material

splint 2.0

Medical resin for hard splints

bg	Инструкция за употреба	5
cs	Návod k použití	8
da	Brugsanvisning	11
de	Gebrauchsanweisung	14
el	Οδηγίες χρήσης	17
en	Instructions for use	20
es	Modo de empleo	23
et	Kasutusjuhend	26
fi	Käyttöohje	29
fr	Mode d'emploi	32
hr	Upute za uporabu	35
hu	Használati útmutató	38
it	Istruzioni per l'uso	41
lt	Naudojimo instrukcija	44
lv	Lietošanas instrukcija	47
nl	Gebruiksaanwijzing	50
no	Bruksanvisning	53
pl	Instrukcja stosowania	56
pt	Manual de instruções	59
ro	Instrucţiuni de utilizare	62
ru	Инструкция по применению	65
sk	Návod na použitie	68
sl	Navodila za uporabo	71
sv	Bruksanvisning	74
tr	Kullanım kılavuzu	77

Ordering information

Freeprint® splint 2.0 385 clear		Freeprint® ibt 385 clear	
500 g	02080	1000 g	04249
1000 g	02076		

Freeprint® crown 385 500 g		Freeprint® crown 385 1000 g	
A1	02372	A1	02376
A2	02378	A2	02415
A3	02417	A3	02446
A3,5	04709	A3,5	04708
B1	02481	B1	02519
B3	02645	B3	02758
C2	02766	C2	02782
D3	02783	D3	02825
BL	02845	BL	02884

Ordering information

Freeprint® denture impact pink-transparent

1000 g	04436
--------	-------

Freeprint® denture impact pink

1000 g	04437
--------	-------

Freeprint® ortho 385 clear

1000 g	03989
--------	-------

Freeprint® tray 2.0 385 green

1000 g	02505
--------	-------

Freeprint® temp 385 500 g

A1	04058
----	-------

A2	04059
----	-------

A3	04060
----	-------

Freeprint® temp 385 1000 g

A1	04062
----	-------

A2	04063
----	-------

A3	04064
----	-------

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Resin for dental 3D printing

Indication

Splints, fixation and transfer keys

Patient Target Group

Persons being treated in the context of a dental procedure.

Intended Users

Dentist, dental technician

Suitable for the following DLP printers/cleaning/post-exposure

see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Minimum material thickness for the design is 1.5 mm
- Polish surface mechanically.
- Processing temperature (at 23 °C ± 2 °C / 73 °F ± 36 °F).

Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- MMA-free.
- Refer to the relevant safety data sheet for hazard and safety information.

Cleaning

By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards. Cleaning tablets can be used following the manufacturer's instructions.

Note

- Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- Do not use any acids or solvents.

Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Please follow the instructions on the safety data sheet!

For users and/or patients

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- Store Freeprint® splint 2.0 in a dry place 15 °C - 28 °C / 59 °F - 83 °F and away from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxide.


Some ingredients of Freeprint® splint 2.0 may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Freeprint® splint 2.0 only insert intraorally in completely polymerised state.

Adverse effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage**Application**

At 23 °C ± 2 °C / 73 °F ± 36 °F

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

① Construction process

Generation of a Print Job complying with machine and material parameters

② Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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detax GmbH