detax

3D Freeprint® Material

temp

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Temporary temps and bridges

Patient target group

Persons being treated in the context of a dental procedure.

Intended users

Dentist, dental technician

Suitable for the following DLP printers/cleaning/post-exposure

see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- Freeprint® temp is applied for the production of temporary c & b in the anterior and posterior area and can be used for single temps and bridges up to 4 elements containing one pontic. The occlusal minimum thickness amounts to 1.5 mm, circularly up to 0.8 mm. The connector cross-section for pontics in the anterior area is 12 mm² minimum, in the posterior area 15 mm². Any non-compliance of these requirements can lead to an undesirable outcome.
- A layer thickness of 50 µm is suitable for the printing process.
- Freeprint® temp can be shaded individually by means of our smartrepair® System.
- Polish the surface mechanically. Prepolish by means of rotating brushes and prepolishing paste, highgloss by means of buffing wheels and highshine polishing paste for resins.
- The temporaries can be cemented by means of usual temporary cements, such as tempolink® (eugenol-free). Temporary cements containing eugenol must be removed completely after the removal of the provisional as otherwise the setting of subsequently used permanent cements could be affected.
- Processing temperature 23 °C ± 2 °C / 73 °F ± 36 °F.

Safety information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Please follow the instructions on the safety data sheet!

For users and/or patients

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- Store Freeprint® temp in a dry place (at 15 °C - 28 °C / 59 °F - 83 °F) and away from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxides.

Some ingredients of Freeprint® temp may cause allergic reactions in predisposed persons.

In such cases refrain from using the product. Freeprint® temp only insert intraorally in completely polymerised state.

Adverse effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

① Construction process

Generation of a Print Job complying with machine and material parameters

② Post-processing


After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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detax GmbH