







# Annex 1






## Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Asiga Max	385 nm	6.1 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
	Asiga Max 2	385 nm	7.0 mW/cm <sup>2</sup>	Composer 2.0.8	Detax_Freeprint tryin 385_5
	Asiga Pico 2	385 nm	20 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
	Asiga PRO 2	385 nm	5.7 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
	Asiga PRO 4K	385 nm	7.0 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
	Asiga Ultra	385 nm	6.6 mW/cm <sup>2</sup>	Composer 2.0.8	Detax_Freeprint tryin 385_5

\*The set of parameters includes all relevant material and printer specific information

# Annex 1



## Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Formlabs Form 4B	405 nm	7.0 mW/cm <sup>2</sup>	PreForm V3.44.0.471	detax Freeprint tryin VO1
	Ivoclar PrograPrint PR5	388 nm	16 mW/cm <sup>2</sup>	PrograPrint CAM 1.1.10.1	Detax FREEPRINT tryin
	Microlay Versus	385 nm	4.3 mW/cm <sup>2</sup>	Microform 1.0.3.7	DETAX Freeprint Tryin 385 100 microns v5.3
	Miicraft Ultra Series	385 nm	5.7 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint tryin 385_100
	Rapidshape D10+/D20+/ D30+/D40+	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint-tryin 385

\*The set of parameters includes all relevant material and printer specific information

# Annex 1

## Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Rapidshape Pro20 (Dental)/ Pro30 (Dental)	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint-tryin 385
	Way2- Production SolFlex Series	385 nm	8.0 mW/cm <sup>2</sup>	Netfabb 2020	Freeprint tryin 385

\*The set of parameters includes all relevant material and printer specific information

# Cleaning Equipment

---

## Cleaning unit

Manufacturer, Model	Cleaning process
---------------------	------------------

---

<b>Ivoclar PrograPrint Clean</b>	<p>Clean the parts with isopropyl alcohol (purity <math>\geq 98\%</math>) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.</p> <p>The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity <math>\geq 98\%</math>) for 3 minutes.</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.</p>
--------------------------------------	--

---

<b>Rapidshape RS wash</b>	<p>Use the following settings: DETAX Freeprint-tryin 385</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues.</p> <p>Then blow off with compressed air.</p>
-------------------------------	--

---

<b>Ultrasonic bath Bandelin Sonorex</b>	<p>Clean the parts with isopropyl alcohol (purity <math>\geq 98\%</math>) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.</p> <p>The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity <math>\geq 98\%</math>) for 3 minutes.</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.</p>
---	--

---



detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany  
+49 7243 510 0 · post@detax.com · detax.com

03/2026

## Light curing Equipment

---

### Light Curing unit

Manufacturer, Model	Curing process
Asiga Cure	Use the following settings: Detax_Freeprint tryin 385_6
Dentalfarm Photopol	2 × 3 min, progressive + N2, turn around components after 3 min
Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 100 %; Duration = 120 s Post curing B: Wavelength = 460 nm; Intensity = 100 %; Duration = 120 s
NK Optik Otoflash G171	2 × 2000 flashes under inert gas, turn around components after 2000 flashes
NK Optik Otoflash 250/500	4000 flashes under inert gas @15 Hz
Rapidshape RS cure	Use the following settings: DETAX Freeprint-tryin 385
Rapidshape RS cure XL	Use the following settings: DETAX Freeprint-tryin 385

---



detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany  
+49 7243 510 0 · post@detax.com · detax.com

03/2026