PERMALINER

Denture Reliner

Definitive Prothesenunterfütterung
Soft denture relining silicone

detax

de	Gebrauchsanweisung	4	no	Bruksanvisning	114	fi	Käyttötarkoitus.....	224
en	Instructions for use	14	pl	Instrukcja stosowania	124	et	Kasutusjuhend.....	234
fr	Mode d'emploi	24	pt	Manual de instruções.....	134	bg	Инструкция за употреба...	244
es	Modo de empleo	34	ro	Instrucțiuni de utilizare	144			
it	Istruzioni per l'uso.....	44	sv	Bruksanvisning	154			
da	Brugsanvisning.....	54	sk	Návod na použitie	164			
el	Οδηγίες χρήσης.....	64	sl	Navodila za uporabo	174			
hr	Upute za uporabu	74	cs	Návod k použití	184			
lv	Lietošanas instrukcija.....	84	hu	Használati útmutató.....	194			
lt	Naudojimo instrukcija	94	ru	Инструкция по применению..	204			
nl	Gebruiksaanwijzing	104	tr	Kullanım Talimatları	214			

Ordering information

PERMALINER

jar of 30 g.....02092

jar of 70 g.....02146

Cutters

Set of 3 pcs.....03101

Molloplast® Pre-Polisher.....03153

10 pcs. pre-polishers

2 pcs. supporting discs

1 pc. mandrel

Grinding sleeves

1 pc. grinding sleeve holder, pointed....02616

10 pcs. grinding sleeves, pointed.....02615

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Intended Use:

Definitive relining of dental prostheses

Patient Target Group:

Persons being treated in the context of a dental procedure.

Intended Users:

Dental technician

Indication:

Suitable for all PMMA-based plastic materials.

- ▶ For fabrication of permanently soft relinings in the dental laboratory.
- ▶ Elimination of pressure points in mandibular and maxillary areas, adhesion problems and flabby ridges.
- ▶ For post-damming the A line, as an elastic margin on maxillary dentures and for securing the position of the dentures by increased adhesion.
- ▶ For elastic covering of sharp-ridged alveolar processes.
- ▶ For covering implant areas and as overdentures.
- ▶ To cover larger defects caused by malformations or operations.

Relining of new dentures with PERMALINER

Plaster cast

- ▶ First, the model is manufactured using type 3 plaster.

Spacer

- ▶ The spacer is ideally produced by deep drawing. Its thickness should be between 1 mm and 3 mm; it is determined by the dentist. The spacer can also be made with cold-curing polymer, e.g. from tray material, and trimmed to the desired thickness.

Investing

- ▶ Invest the wax pattern with spacer in a flask for the press technique. Insulation of the teeth in the counter with **FLEXISTONE®** or **FLEXISTONE® Plus** is recommended. Place the flask in boiling water for approx. 4 minutes, then open and remove all traces of wax carefully. Remove spacer first.

Insulation

- ▶ Apply conventional alginate insulation. Then return the spacer to its position.

Pressing and pre-polymerization of the acrylic

- ▶ The acrylic is mixed and packed according to the instructions for use. A polyethylene foil is placed between the acrylic and the spacer for the pressing process. Pre-press the flask at 1.5 tons at intervals of approx. 10–15 minutes. After opening the flask, remove the polyethylene foil and the excess acrylic. Close the flask again with the spacer and the polyethylene foil and press it once more with the clamp. Close the clamp, place the flask in cold water and bring to the boil. Boil for about 30 minutes, then let it cool down. This pre-polymerization avoids possible reactions with the acrylic monomer. Remove the flask from the clamp and open it. Remove spacer and polyethylene foil. Before the packing of **PERMALINER**, the acrylic must have reached a high degree of stiffness to avoid deformation.

Packing of PERMALINER

- ▶ Withdraw the ready-to-use **PERMALINER** from the jar with a clean spatula, place it on the acrylic and press it on slightly. Place a new polyethylene foil between **PERMALINER** and the flask counter. Close the flask and pre-press at intervals with approx. 1.5 tons. After opening the flask, remove polyethylene foil and **PERMALINER** excess. Close the flask again and perform the final pressing with approx. 3 tons. For polymerization, place the flask in the clamp under pressure.

Polymerization

- ▶ Place the clamped flask in a cold water bath and heat slowly up to 100 °C / 212 °F. Polymerization time: **Boil for 1 h at 100 °C / 212 °F.** Then allow the flask to cool down slowly. Do not quench with cold water.

Polymerization in microwave device

- ▶ Polymerization takes place in suitable commercially available microwave flasks (plastic, suitable for pressing) at 850 watts within 10 minutes. Microwave-compatible dental plasters must be used for investing. Investing and polymerization must be carried out within 8 hours. Units without a rotary plate may necessitate gradual polymerization on all sides by turning of the flask. De-flask after self-cooling down to room temperature.

Trimming of the denture

- ▶ Trim the acrylic as usual. **PERMALINER** can be machined with milling cutters or grinding sleeves (15,000 – 20,000 RPM) without heat build-up. Rough spots can be removed with polishing discs.

Polishing of the acrylic

- ▶ Pre-polish and mirror-polish acrylic in the conventional way. **PERMALINER** cannot be mirror-polished. **Lustrol** gloss varnish is recommended as an alternative. The dosage is 1:1. Before application, the base and

functional borders must be clean and dry. Apply **Lustrol** using a brush and let it dry for 15–20 minutes. **Lustrol** is suitable for use in combination with **Molloplast® B**, **PERMALINER**, **mollosil®**.

Relining of an existing acrylic denture with PERMALINER

- ▶ In the case of an existing denture that is to be constantly softly relined with **PERMALINER** the dentist makes a new impression of the base with the existing denture. Medium- or light-bodied impression materials are used for producing the impression. Model fabrication, investing and insulation are carried out as described above.

The acrylic base

- ▶ Take the denture from the counter and remove the impression material from the denture. Reduce the base with milling by the desired **PERMALINER** thickness (1.0–3.0 mm). For a seamless transition from acrylic to **PERMALINER** it is recommendable to cut an even circular shoulder in the transition area into the acrylic in the desired **PERMALINER** thickness. Completely remove any existing soft relining material.
- ▶ **Packing, pressing, polymerizing** and trimming of **PERMALINER** are carried out as described above.

Recommendations to avoid incorrect processing

1. Blistering:

- ▶ Insufficient pressure during pressing.
- ▶ Flask placed in hot – rather than cold – water bath.
- ▶ Packing of **PERMALINER** before the acrylic is stiff enough.
- ▶ Wetting of the acrylic with monomer liquid before application of **PERMALINER**.

2. Insufficient or no bonding:

- ▶ Prosthetic base plate not roughened deeply enough.
- ▶ The denture was wetted with alginate insulation.
- ▶ Hands and denture were not grease-free.
- ▶ Excessive heat build-up during trimming can cause **PERMALINER** to peel off at the borders.

Recommendations for repairing **PERMALINER**

- ▶ **PERMALINER** can be repaired and supplemented with itself. Invest the denture for this purpose. Roughen existing **PERMALINER** and plastic material well. For supplementation in the repair area, place new **PERMALINER** in the flask and press (trial pressing with foil possible). Then polymerize according to instructions.

Safety Information

- ▶ Only for the specified intended use by trained specialists.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

Notes

- ▶ Detax shall not be held liable for any damage caused by misuse.
- ▶ Read and understand the safety data sheet!

For users and / or patients

Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the country in which the user and/or patient resides.

Contraindication

Contains organic peroxides.

Ingredients of **PERMALINER** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply **PERMALINER** intraorally only when fully polymerized.

Side Effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Symbols

Cleaning recommendations for your permanently softly relined denture

For the patient

Please make a copy and give it to the patient

Your new **PERMALINER** relining has many benefits. It looks like natural gingiva, is comfortable to wear, and is resistant to mechanical pressure. To ensure that you can enjoy your **PERMALINER** relinings for a long time, you should observe the following care instructions:

- ▶ Remove the denture after each meal and clean it thoroughly with water.
- ▶ Especially after drinking fruit or vegetable juices, rinse the mouth well with water.
- ▶ 1 × a day the denture must be cleaned thoroughly, removing all plaque. Liquid soap and a soft toothbrush are suitable for this.
- ▶ Commercially available denture cleaning agents may also be used. However, in order to avoid bleaching, the denture should not remain in the solution for longer than 10–15 minutes.
- ▶ At night the relined denture should be removed, cleaned, and stored in a dry place.
- ▶ Frequent consumption of tea, coffee or wine can cause discoloration.
Special note to smokers: Tobacco discolors the relining irreversibly!

Additional care with Detax “good morning spray”

Detax good morning spray improves oral hygiene. Colonization with oral-specific bacteria is delayed. Wet the cleaned underside of the denture in the morning with 2 – 3 spray strokes before inserting it, and insert it immediately.

Note: Does not replace the cleaning of the denture.

Also clean your teeth so as to remove all plaque according to the instructions of your dentist. You benefit in any case.

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detax GmbH