Esthetic Mask

DE

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Esthetic Mask

GB
Instructions for use

Flexible gingival mask for the fabrication of master models in the dental technique, 1:1 A-silicone-based, light bodied, cold curing

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Fig. 6


Fig. 7

Fig. 8

1. Preparation of the silicone mold (indirect)

Before sawing the model, build a mold by creating a bank around all parts of the model to be masked, using addition or condensation curing putty (e.g. **blue eco** or **compact lab putty**) (Fig. 1). An existing impression, if made of silicone, can also be used as a bank after separating it from its tray. Be sure to remove any undesired protrusions from areas outside of the mask area on the bank to later allow better realignment on the model.

2. Preparation of the model

Grind away a layer of material in those areas to be modelled in **Esthetic Mask** (Fig. 2). Be sure to remove enough material so as to allow sufficient mask thickness. Stumps can be sawed and worked on either before or after mask preparation; if done before mask preparation, be sure to block-out saw lines.

3. Preparing silicone mold for assembly

Cut several injection holes into the silicone mold (from the palatal resp. lingual side). Then, drill air vents into the mold with a rose-head burr, placing them at the highest points of the desired mask (Fig. 3). Clean silicone mold and coat **with separating liquid**. Let dry for 30 seconds, then remount the mold onto the model, check to ensure satisfactory fit, and fasten.

4. Fabrication by direct method

Mark the edge of the gingival areas for the impression with soft wax or similar. Then insulate this area by spraying lightly with silicone Separating Agent. Drying time 30 sec. **Separating Agent** (separating liquid) is suitable for A- and C-silicones and polyether impression material. Please note by the application of **Esthetic Mask**: Inject directly around the impression or transfer posts and leave the tip of the mixer in the material during dispensing (Fig. 7). Dispense with steady pressure. Ensure that the **Esthetic Mask** has the desired dimensions. Continue model fabrication only after the **Esthetic Mask** has set completely. Any separating medium residues on the laboratory posts can be removed with isopropyl alcohol. Insulation of the mask against plaster or resins is not necessary.

5. Mixing and casting

Extrude equal lengths of base material and catalyst onto a glass plate or mixing pad (Fig. 4) and quickly homogenise with a spatula (30 sec.). An even, bubble-free mixture is best attained by repeatedly spreading the material out on the mixing pad and scraping it off again with the spatula. Holding the disposable syringe at an angle, fill with material while jiggling (Fig. 5), insert plunger at back of syringe and quickly fill mold (Fig. 6). Inspection: air vent channels must be partially filled with silicone. **Esthetic Mask** will set in 8 - 10 minutes and can then be removed from the mold (Fig. 8).

6. Trimming and finishing

When the mask has been fabricated on the model, remove the mould carefully from the model, trim any casting residues from the air vents with a scalpel and remove the mask. Remove any flash with a scalpel or sharp scissors. Further processing of the removed mask is possible with suitable rotating tools (Molloplast® cutters, grinding sleeves with 15,000 - 25,000 r.p.m., tungsten carbide drills and arbors).

Important hints

- Do not interchange tube closures.
- Make sure that only base and catalyst components from the same lot number are used.
- Syringes can be cleaned and re-used after **Esthetic Mask** has set.
- **Esthetic Mask** can also be pressure-cured at 2.5 bar (dry pressure).
- **Esthetic Mask** is not subject to warping or other dimensional changes.
- **Esthetic Mask** can be cleaned in a standard steam-cleaning apparatus.
- Contact with latex examination gloves or other latex material may affect the setting reaction of **Esthetic Mask**.
- For demo and presentation models **Lustrol** gloss varnish can be applied for a high-gloss surface.
- Please follow the instructions of the safety data sheets!

Only for dental use by qualified personnel.

approx. 2 min.

Working time

approx. 8 - 10 min.

Setting time

Indications for use:

- Preparation of gingival masks on master models and saw models
- Crown and bridge technique
- Colour matching for veneering techniques
- Combination crown techniques
- Tapered and telescoping crown technique
- Implant construction

Technical data:

- **Mixed volume:** 2 x 50 ml (tubes)
- **Mixing ratio:** 1:1
- **Colour code:** base: reddish brown catalyst: white
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 2 min.*
- **Setting time:** approx. 8 - 10 min.*
- **Recovery from deformation:** approx. 99%
- **Linear dimensional change:** approx. 0.3%
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from initiation of mixing a 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

Esthetic Mask
Standard packing 02340

Base, tube of 50 ml
Catalyst, tube of 50 ml
Separating liquid for silicones, vaporizer of 15 ml
2 Disposable syringes

Separating Agent 02690
vaporizer of 15 ml

Grinding sleeves 02615
pointed, 10 units

Grinding sleeve holder 02616
pointed, 1 unit

Esthetic Mask automix
Standard packing 03113
cartridges of 2 x 50 ml
10 mixing cannulas, yellow
15 ml Separating Agent

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