Esthetic Mask automix

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Esthetic Mask automix

EN
Instructions for use

Flexible gingival mask for dental laboratory implant models, A-silicone based 1:1, light bodied, normal setting


Fig. 1

Fig. 2

Fig. 3

Fig. 4


Fig. 5

Fig. 6


Fig. 7

Fig. 8

1. Preparation

a) Fabrication by direct method, application to impression

Mark the edge of the gingival areas for the impression with soft wax or similar. Then insulate this area by spraying lightly with silicone **Separating Agent**. Drying time 30 sec. **Separating Agent** (separating liquid) is suitable for A- and C-silicones and poly-ether impression material. Please note by the application of **Esthetic Mask automix**: Inject directly around the impression or transfer posts and leave the tip of the mixer in the material during dispensing (Fig. 1). Dispense with steady pressure. Ensure that the **Esthetic Mask automix** has the desired dimensions. Continue model fabrication only after the **Esthetic Mask automix** has set completely. Any separating medium residues on the laboratory posts can be removed with isopropyl alcohol. Insulation of the mask against plaster or resins is not necessary.

b) Fabrication by indirect method, on the model

To fabricate the **Esthetic Mask automix** first make a silicone mould of addition- or condensation curing putty (e.g. **blue eco** or **compact lab putty**), which covers all of the areas of the unsawn model to be reproduced (Fig. 2). An available impression can be used as mould after removing the tray. Remove any areas in the index or impression that are in the way outside the gingival mask region so that it can be repositioned easily on the model. Reduce the regions of the model that are to be reproduced by **Esthetic Mask automix** by milling (Fig. 3). Take care to remove enough material to ensure the later thickness of the mask. Sawing and working the dies can take place before or after mask fabrication depending on the situation. If sawing is performed beforehand, block out the saw cuts. Drill one or more injection openings (palatal/lingual) in the silicone. Drill air vent channels through the mould with a round bur starting from the highest points of the mask (Fig. 4), clean the mould and insulate it on the inside by spraying lightly with silicone **Separating Agent**. Drying time 30 sec. Put the mould back on the model, check that it is seated exactly and fix it. Inject material quickly through the injection opening of the mould with the cartridge (Fig. 5). Inspection: air vent channels must be partially filled with silicone. Demould the mask only after complete setting.

2. Mixing and applying

The cartridge is extruded with the Automix2 system gun (Fig. 6). After inserting the cartridge into the mixing gun, remove the cartridge cap by twisting off. Express a small amount of material until silicone is dispensed evenly from the two openings (Fig. 7). To attach the mixing cannula note the guides on the cannula and cartridge. Rotate anti-clockwise to a stop (Fig. 8). Dosing can now take place individually. Dispense the material with steady pressure. After use leave the mixing cannula on the cartridge until the next use.

3. Processing

When the mask has been fabricated on the model, remove the mould carefully from the model, trim any casting residues from the air vents with a scalpel and remove the mask. Remove any flash with a scalpel or sharp scissors. Further processing of the removed mask is possible with suitable rotating tools (Molloplast® cutters, grinding sleeves with 15,000 - 25,000 r.p.m., tungsten carbide drills and arbors).

Important working hints

- **Esthetic Mask automix** can also be pressure-cured at 2.5 bar (dry pressure).
- **Esthetic Mask automix** is not subject to warping or other dimensional changes.
- **Esthetic Mask automix** can be cleaned in a standard steam-cleaning apparatus.
- Contact with latex examination gloves or other latex material may affect the setting reaction of **Esthetic Mask automix**.
- For demo and presentation models **Lustrol** gloss varnish can be applied for a high-gloss surface.
- Please follow the instructions of the safety data sheets!
- Only for dental use by qualified personnel.

approx. 1 min. 30 sec.

Working time

approx. 7 min.

Setting time

Indications for use:

- Gingival masks on master models and saw models
- C & B technique
- Natural colour matching for veneering techniques
- Combination crown techniques
- Tapered & telescoping crown technique
- Implant construction

Technical data:

- **Mixing volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: reddish brown catalyst: white
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 1 min. 30 sec.*
- **Setting time:** approx. 7 min.*
- **Recovery from deformation:** approx. 99,5 %
- **Linear dimensional change:** approx. 0,2 %
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

Esthetic Mask automix
Standard packing 03113
cartridges of 2 x 50 ml
10 mixing cannulas, yellow
15 ml Separating Agent

Mixing gun 02699
Automix2, 1 pc.

Mixing cannulas 02706
yellow, 48 pcs.

Intra-Oral-Tips 02712
yellow, 96 pcs.

Esthetic Mask
Standard packing 02340
base, tube of 50 ml
catalyst, tube of 50 ml
15 ml Separating Agent
2 disposable syringes

Implant Mask
Standard packing 02529
cartridges of 2 x 10 ml
6 mixing cannulas,
light green, mini-mix
15 ml Separating Agent

Separating Agent 02690
separating liquid for sili-
cones, vaporizer of 15 ml

blue eco scan
Standard packing 03114
base, jar of 700 g
catalyst, jar of 700 g
2 measuring scoops

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