
Molloplast® B

permanently soft

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Instructions for use

Intended Use:

Definitive relining of dental prostheses

Patient Target Group:

Persons being treated in the context of a dental procedure.

Intended Users:

Dental technician

Product description:

Heat-curing one-component silicone, permanently soft relining material, suitable for denture plastics based on (meth)acrylate.

Indications:

Sharp-ridged alveolar processes with thin tegument. Cushioning of the denture against the sharply defined linea mylohyoidea, against the foramen mentale or other prominent points of the alveolar process, such as on the trigonum retromolare or the tubera. Cushioning against the sutura mediana. For post-damming of the A line and stabilizing of the position of the maxillary and mandibular dentures by increased adhesion. To cover larger defects caused by malformations or surgeries with obturators.

Relining with Molloplast® B with simultaneous fabrication of a new denture

Plaster cast

First pour a master cast from dental stone class III (Fig. 1).

The spacer

is ideally produced by deep drawing. Its thickness should be between 1 mm and 3 mm and is determined by the dentist. If no thermoforming unit is available in the laboratory, the spacer can be made with cold-curing polymer and later trimmed to the desired thickness (Fig. 2). Set up teeth and try-in with spacer in the patient's mouth (Fig. 3).

Investing

After the try-in, invest the wax pattern with spacer into the flask with counter. Insulation of the teeth with FLEXISTONE® or FLEXISTONE® Plus is recommended. Place the flask in boiling water for approx. 4 minutes, then open it and remove all traces of wax.

Important:

Wax residues must not penetrate or adhere to the plaster.

Insulation

Apply an alginate-based insulating liquid. Then reposition the spacer.

Mixing of the acrylic

Mix the acrylic according to the instructions for use and pack it into the hand-warm flask. Place a polyethylene

foil between the acrylic and the spacer for the pressing process (Fig. 4)

Pressing

Pre-press the flask at intervals of approx. 10 – 15 min at approx. 100 kp. After opening the flask, remove the polyethylene foil and the excess acrylic. Then close the flask with the spacer and the polyethylene foil again and press it once more with the clamp. Close the clamp, place the flask in cold water and bring to the boil. Boil for about 30 minutes, then let it cool down. This pre-polymerization avoids possible reactions with the acrylic monomer.

Packing of Molloplast® B

Remove the flask from the clamp and open it. Remove spacer and polyethylene foil.

Important:

In order to avoid deformation, the acrylic must have reached a high degree of stiffness before the packing of Molloplast® B (Fig. 5). Withdraw the ready-to-use Molloplast® B from the jar with a clean plastic spatula and pack it onto the pre-pressed acrylic. Place a new polyethylene foil between Molloplast® B and the flask counter. Close the flask and pre-press at intervals with approx. 100 kp. After opening the flask, remove the polyethylene foil and Molloplast® B excess. Close the flask again and press for approx. 10 – 15 minutes at 100 kp. For polymerization, place the flask into the clamp.

Polymerization

Place the flask in a cold water bath and slowly heat the water to 100 °C / 212 °F.

Polymerization time: 2 h at 100 °C / 212 °F.

Let the flask cool down slowly. Do not quench with cold water!

Molloplast® B polymerization in the microwave device

Polymerization is carried out in suitable commercially available microwave flasks at 850 watts within 10 minutes. Microwave-compatible dental plaster must be used for investing. Investing and polymerization must be carried out within 8 hours. Units without a rotary plate may necessitate gradual polymerization on all sides by turning of the flask. De-flask after self-cooling down to room temperature.

Trimming of the denture

Trim the acrylic as usual. The prosthetic base plate and functional borders made of Molloplast® B can be trimmed with Molloplast® special cutters or grinding

sleeves (15,000 – 20,000 RPM) easily and without heat build-up. Use Molloplast® pre-polishing discs for smoothing of rough spots (Fig. 6).

Polishing of the acrylic

Pre-polish the acrylic in the conventional way with pumice powder, brush or buffering tool, and polish to a high gloss using commercially available materials. Mechanical high-gloss polishing is not possible with Molloplast® B; to this end, Lustrol is recommended. The dosage is 1:1. Before applying the liquid polish, base and functional borders must be clean and dry. Apply Lustrol with a brush and allow it to dry for approx. 10 – 15 minutes. The high-gloss polishing is finished (Fig. 7). The newly fabricated and permanently softly relined denture is now handed over to the dentist for insertion.

Initial relining of an existing acrylic denture with Molloplast® B

In the case of an existing denture that is to be permanently softly relined with Molloplast® B, the dentist makes a new impression of the base with the existing denture. The impression is taken with a medium- or light-bodied impression material (Fig. 8).

Plaster cast

A master cast of the new impression is made with dental stone class III. After the curing of the plaster model, it is invested into the flask with counter. It is recommended to cover the teeth and a part of the denture surface with silicone (FLEXISTONE® or FLEXISTONE® Plus) before

creation of the flask counter, as this makes it easier to remove the denture from the plaster counter. After the curing of the plaster, open the flask and remove the impression material from the denture.

The acrylic base

is reduced by approx. 1 – 3 mm with a grinding stone.

Important:

For a seamless transition from acrylic to Molloplast® B it is recommendable to grind grooves evenly into the functional border using a wheel or fissure bur. Then use a brush to apply Primo adhesive carefully once or twice onto the entire acrylic base. Let Primo air-dry for 60 – 90 minutes before applying Molloplast® B (Fig. 9). Insulate the plaster model with a commercially available alginate-based insulating liquid.

Packing of Molloplast® B

takes place with a hand-warm flask. Remove the ready-to-use Molloplast® B from the jar with a clean plastic spatula and apply it evenly to the prosthetic base plate. Cover with a polyethylene foil, close the flask and place it under the press.

Pressing

Press for approx. 4 min, open flask, remove polyethylene foil and excess Molloplast® B. Check whether enough Molloplast® B has been applied to the base. Apply more if necessary. Close the flask again and leave it under the press for 15 minutes at a pressure of approx. 100 kp. Place flask in the clamp and tighten.

Polymerization

Place the flask in a cold water bath and slowly heat the water to 100 °C / 212 °F.

Polymerization time: 2 h at 100 °C / 212 °F.

Let the flask cool down slowly. Do not quench with cold water!

Molloplast® B polymerization in the microwave device

Polymerization is carried out in suitable commercially available microwave flasks at 850 watts within 10 minutes. Microwave-compatible dental plaster must be used for investing. Investing and polymerization must be carried out within 8 hours. Units without a rotary plate may necessitate gradual polymerization on all sides by turning of the flask. De-flask after self-cooling down to room temperature.

Trimming of the denture

Trim the acrylic as usual. The prosthetic base plate and functional borders can be trimmed with Molloplast® special cutters or grinding sleeves (15,000 – 20,000 RPM) easily and without heat build-up. Use Molloplast® pre-polishing discs for smoothing of rough spots (Fig. 6).

Polishing of the acrylic

Pre-polish the acrylic in the conventional way with pumice powder, brush or buffering tool, and polish to a high gloss using commercially available materials. Mechanical high-gloss polishing is not possible with Molloplast® B; to this end, Lustrol is recommended. The dosage is 1:1. Before applying the liquid polish, base and functional borders must be clean and dry. Apply Lustrol with a brush and allow it to dry for approx. 10 – 15 minutes. The high-gloss polishing is finished (Fig. 7).

Primo adhesive

The chemical bond between acrylic and Molloplast® B is established by the Primo adhesive. This applies only to denture material already polymerized. Prior to the application of Primo adhesive, the prosthetic base plate must be prepared according to Molloplast® B instructions for use. Primo adhesive is applied to the cleaned and dried prosthetic base plate once or twice over the entire surface. Note: The entire contact surface acrylic – Molloplast® B must be covered. After a drying time of approx. 60 – 90 minutes, Molloplast® B is pressed on **and then polymerized according to instructions.**

Recommendations to avoid incorrect processing

1. Air inclusions in polymerized dentures:

- a) Not enough pressure during pressing.
- b) Flask placed in hot – rather than cold – water bath.
- c) Packing of Molloplast® B before acrylic is stiff.
- d) Wetting of the acrylic with monomer liquid before application of Molloplast® B.

2. Insufficient or no bonding of Molloplast® B to the acrylic:

- a) Too few retentions have been cut into the existing denture base, or Primo adhesive has not been applied.
- b) The acrylic denture has been wetted with alginate-based insulating liquid.
- c) Hands and working material were not free of grease.
- d) Excessive heat build-up during trimming can cause Molloplast® B to peel off at the functional borders. Use only Molloplast® special cutters and Molloplast® pre-polishers for trimming.

Recommendations for repair of Molloplast® B

Molloplast® B can be repaired with Molloplast® B in case of breakage of a denture or for partial supplementation. Invest the denture as usual. Grind off approx. 1 mm of old material in the area to be repaired (using Molloplast® cutters). In case of a fracture to be repaired, open the area in a V shape, starting from the base. Apply new Molloplast® B for supplementation to the repair area of the denture invested in the flask and put it under the press. (Test pressing with foil is possible.) Then polymerize according to instructions.

Do not use Primo adhesive!

Safety Information

- Only for the specified intended use by trained specialists.
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- In case of accidental contact, rinse eyes immediately with running water for at least 15 minutes, consult an ophthalmologist.
- Primo adhesive: Do not inhale vapors. Keep away from sources of ignition.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- DETAX shall not be held liable for any damage caused by misuse.

- Read and understand the safety data sheet!

For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.de and to the competent public authority of the country in which the user and/or patient resides.

Contraindication

Contains organic peroxides.

Ingredients of Molloplast® B may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply Molloplast® B intraorally only when it is fully polymerized.

Side Effects


Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Symbols

MD Medical device

see base

Cleaning recommendations for your permanently softly relined denture

**For the patient
Please make a copy and give it to
the patient**

Your new Molloplast® B relining has many benefits. It looks like natural gingiva, is comfortable to wear, and is resistant to mechanical pressure.

To ensure that you can enjoy your Molloplast® B relinings for a long time, you should observe the following care instructions:

- Remove the denture after each meal and clean it thoroughly with water.
- Especially after drinking fruit or vegetable juices, rinse the mouth well with water.

- 1 x a day the denture must be cleaned thoroughly, removing all plaque. Liquid soap and a soft toothbrush are suitable for this.
- Commercially available denture cleaning agents may also be used. However, in order to avoid bleaching, the denture should not remain in the solution for longer than 10–15 minutes.
- At night the relined denture should be removed, cleaned, and stored in a dry place.

Frequent consumption of tea, coffee or wine can cause discoloration.
Special note to smokers: Tobacco discolors the relining irreversibly!

Also clean your teeth so as to remove all plaque according to the instructions of your dentist. You benefit in any case.

Additional care with DETAX "good morning spray"

DETAX good morning spray improves oral hygiene. Colonization with oral-specific bacteria is delayed. Wet the cleaned underside of the denture in the morning with 2 – 3 spray strokes before inserting it, and insert it immediately.

Note: Does not replace the cleaning of the denture.

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



Ordering information	Lustrol	FLEXISTONE®	FLEXISTONE® Plus
Molloplast® B <u>Standard packing</u> jar of 45 g	03000 <u>Standard packing</u> bottle of 6 ml gloss varnish bottle of 6 ml catalyst 2 pipettes	03008 <u>Standard packing</u> 190 ml bottle (350 g) 30 ml catalyst, bottle 1 mixing bowl	03009 <u>Standard packing</u> 2 x 160 ml bottles base + catalyst 1 mixing bowl
<u>Large packing</u> jar of 170 g	03001 <u>Molloplast® -Cutters</u> 3 pcs. set 2 pcs. set	03101 <u>Large packing</u> 5 x 190 ml bottles (350 g) 03120 <u>Molloplast® -pre-polishers</u> 10 pcs. + 1 mandres	03010 <u>Eco packing</u> 2 x 1000 ml bottles base + catalyst 03153 <u>Catalyst</u> 30 ml, bottle
<u>Combi packing</u> jar of 45 g Dose bottle of 15 ml Primo Adhesive	03002 <u>Grinding sleeves</u> 1 St. Grinding sleeve holder pointed	 02616 10 St. Grinding sleeves pointed	 03165 <u>good morning spray</u> Standard packing 100 ml, spray bottle
Primo Adhesive bottle of 15 ml with brush	03004 <u>Grinding sleeves</u> 1 St. Grinding sleeve holder pointed	 02615	 <u>Eco packing</u> 5 x 100 ml, spray bottle
Molloplast® Remover 10 ml dropping bottle	02206		02811

Molloplast® B

permanently soft

9

Application Clip on YouTube!

DIN EN ISO 10139 - 2, Typ A

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).