

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Aidite Cameo CPD-100	405 nm	4.3 mW/cm ²	Aidite Dental	EZPRINT Crown Bridge
	Asiga Max	385 nm	6.1 mW/cm ²	Composer 1.2.11	Detax_Freeprint crown_5
	Asiga Max 2	385 nm	7.0 mW/cm ²	Composer 2.0.8	Detax_Freeprint crown_5
	Asiga PRO 4K	385 nm	7.0 mW/cm ²	Composer 1.2.11	Detax_Freeprint crown_5
	Asiga Ultra	385 nm	6.6 mW/cm ²	Composer 2.0.8	Detax_Freeprint crown_5

*The set of parameters includes all relevant material and printer specific information

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Heygears UltraCraft A2D HD	385 nm	1.75 mW/cm ²	Heygears Cloud V4.3.7.6	Detax Crown
	Formlabs Form 4B	405 nm	7.0 mW/cm ²	PreForm V3.44.0.471	detax Freeprint crown V01
	Ivoclar PrograPrint PR5	388 nm	16 mW/cm ²	PrograPrint CAM11.10.1	Detax FREEPRINT crown
	LuxCreo FastPrint.io	385 nm	6 mW/cm ²	LuxFlow v31.0.301229	DETAX freeprint-crown
	Microlay Versus	385 nm	4.3 mW/cm	Microform 1.0.3.7	DETAX Freeprint Crown 385 50 microns v5.6

*The set of parameters includes all relevant material and printer specific information

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Rapidshape ONE	385 nm	3.5 mW/cm ²	Print Studio 1.2.	DETAX Freeprint-crown
	Rapidshape D10+/D20+/ D30+/D40+	385 nm	2.0 mW/cm ²	Netfabb 2020	DETAX Freeprint-crown
	Rapidshape PRO20 (Den- tal)/PRO30 (Dental)	385 nm	2.0 mW/cm ²	Netfabb 2020	DETAX Freeprint-crown
	Rayshape Edge E2	405 nm	5.0 mW/cm ²	ShapePanel	DETAX FREEPRINT crown
	Rayshape Edge Mini	405 nm	5.0 mW/cm ²	ShapeWare 2.0 V0.5.5	DETAX FREEPRINT crown

*The set of parameters includes all relevant material and printer specific information

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Shining 3D AccuFab CEL	405 nm	4.0 mW/cm ²	AccuWare V3.2.1.109	DETAX FREEPRINT crown -0.05
	Shining 3D AccuFab L4D	405 nm	2.0 mW/cm ²	AccuWare V3.2.1.108	DETAX FREEPRINT crown -0.05
	W2P SolFlex 250 UHD	385 nm	5.0 mW/cm ²	SolFab	Detax FREEPRINT crown 385

*The set of parameters includes all relevant material and printer specific information

Cleaning Equipment

Cleaning unit

Manufacturer, Model	Cleaning process
---------------------	------------------

Ivoclar PrograPrint Clean	<p>Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.</p> <p>The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 2 minutes.</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.</p>
--------------------------------------	--

Rapidshape RS wash	<p>Use the following settings: DETAX Freeprint-crown</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues.</p> <p>Then blow off with compressed air.</p>
-------------------------------	--

Ultrasonic bath Bandelin Sonorex	<p>Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 1 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.</p> <p>The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 1 minutes.</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.</p>
---	--



detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany
+49 7243 510 0 · post@detax.com · detax.com

02/2026

Light curing Equipment

Light Curing unit

Manufacturer, Model Curing process

Asiga Cure Use the following settings: Detax_Freeprint crown_6

Dreve PCU LED N2 Intensity: 40%, 6:40 min, vacuum, nitrogen

Heygears PCU 3.0 Use the following settings: Detax Crown UV

Ivoclar PrograPrint Cure Post curing A: Wavelength = 405 nm; Intensity = 45 %; Duration = 5 min
Post curing B: Wavelength = 460 nm; Intensity = 45 %; Duration = 5 min

LuxCreo ILCD Program „PO2“

NK Optik Otoflash G171 2 × 2000 flashes under inert gas, turnaround components after 2000 flashes

NK Optik Otoflash 250/500 4000 flashes under inert gas @15 Hz

Rapidshape RS cure Use the following settings: DETAX Freeprint-crown

Rapidshape RS cure XL Use the following settings: DETAX Freeprint-crown



detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany
+49 7243 510 0 · post@detax.com · detax.com

02/2026