

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufacturer, Model	Light source	Pixel size x,y	Light intensity	Operation Software	Parameter data set*	Curing time	Slicing	Recommended orientation angle (degree)
	Aidite Cameo CPD-100	405 nm	±30 µm	4.3 mW/cm ²	Aidite Dental	EZPRINT Crown Bridge	2.4 s	50 µm	0-90
	Asiga Max	385 nm	±62 µm	6.1 mW/cm ²	Composer 1.2.11	Detax_Freeprint crown_5	2.15 s	50 µm	0-90
	Asiga Max 2	385 nm	±62 µm	7.0 mW/cm ²	Composer 2.0.8	Detax_Freeprint crown_5	2.0 s	50 µm	0-90
	Asiga PRO 4K	385 nm	±56 µm	7.0 mW/cm ²	Composer 1.2.11	Detax_Freeprint crown_5	2.0 s	50 µm	0-90
	Asiga Ultra	385 nm	±50 µm	6.6 mW/cm ²	Composer 2.0.8	Detax_Freeprint crown_5	2.1 s	50 µm	0-90
	Formlabs Form 4B	405 nm	±50 µm	7.0 mW/cm ²	PreForm V3.44.0.471	detax Freeprint crown V01	1.5 s	50 µm	0-90
	Heygears UltraCraft A2D HD	385 nm	±42.5 µm	1.75 mW/cm ²	Heygears Cloud V4.3.7.6	Detax Crown	3.0 s	50 µm	0-90
	Ivoclar PrograPrint PR5	388 nm	±25 µm	16 mW/cm ²	PrograPrint CAM 1.1.10.1	Detax FREEPRINT crown	0.9 s	50 µm	0-90
	LuxCreo Fast-Print.io	385 nm	±40 µm	6.0 mW/cm ²	LuxFlow v31.0.301229	DETAX freeprint-crown	2.15 s	50 µm	0-90
	Microlay Versus	385 nm	±65 µm	4.3 mW/cm ²	Microform 1.0.3.7	DETAX Freeprint Crown 385 50 micron v5.6	1.99 s	50 µm	0-90
	Rapidshape ONE	385 nm	±34 µm	3.5 mW/cm ²	Print Studio 1.2.	DETAX Freeprint-crown	1.8 s	50 µm	0-90
	Rapidshape D10+/D20+/D30+/D40+	385 nm	±34 µm	2.0 mW/cm ²	Netfab 2020	DETAX freeprint-crown	3.6 s	50 µm	0-90
	Rapidshape PRO20 (Dental)/PRO30 (Dental)	385 nm	±34 µm	2.0 mW/cm ²	Netfab 2020	DETAX Freeprint-crown	3.6 s	50 µm	0-90
	Rayshape Edge E2	405 nm	±34.4 µm	5.0 mW/cm ²	ShapePanel	DETAX FREEPRINT crown	1.8 s	50 µm	0-90
	Rayshape Edge Mini	405 nm	±35 µm	5.0 mW/cm ²	ShapeWare 2.0 V0.5.5	DETAX FREEPRINT crown	2.5 s	50 µm	0-90
	Shining 3D AccuFab CEL	405 nm	±35 µm	4.0 mW/cm ²	Accu Ware V3.2.1.109	DETAX FREEPRINT crown -0.05	2.0 s	50 µm	0-90
	Shining 3D AccuFab L4D	405 nm	±50 µm	2.0 mW/cm ²	Accu Ware V3.2.1.108	DETAX FREEPRINT crown -0.05	3.7 s	50 µm	0-90
	W2P SolFlex 250 UHD	385 nm	±50 µm	5.0 mW/cm ²	SolFab	Detax FREEPRINT crown 385	2.4 s	50 µm	0-90

*The set of parameters includes all relevant material- and printer specific information

Cleaning Equipment

Cleaning unit Manufacturer, Model	Cleaning process
Ivoclar PrograPrint Clean	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 2 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Rapidshape RS wash	Use the following settings: DETAX Freeprint-crown Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 1 minute. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 1 minute. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

Light curing Equipment

Light Curing unit Manufacturer, Model	Curing process
Asiga Cure	Use the following settings: Detax_Freeprint crown_6
Dreve PCU LED N2	Intensity: 40%, 6:40 min, vacuum, nitrogen
Heygears PCU 3.0	Use the following settings: Detax Crown UV
Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 45 %; Duration = 5 min Post curing B: Wavelength = 460 nm; Intensity = 45 %; Duration = 5 min
LuxCreo ILCD	Program „PO2“
NK Optik Otofash G171	2 x 2000 flashes under inert gas, turn around components after 2000 flashes
NK Optik Otofash 250/500	4000 flashes under inert gas @15 Hz
Rapidshape RS cure	Use the following settings: DETAX Freeprint-crown
Rapidshape RS cure XL	Use the following settings: DETAX Freeprint-crown

