detax

3D **Freeprint®** Material

ortho

Medical resin for surgical guides

bg	Инструкция за употреба	3
cs	Návod k použití	7
da	Brugsanvisning	11
de	Gebrauchsanweisung	15
el	Οδηγίες χρήσης	19
en	Instructions for use	23
es	Modo de empleo	27
et	Kasutusjuhend	31
fi	Käyttöohje	35
fr	Mode d'emploi	39
hr	Upute za uporabu	43
hu	Használati útmutató	47
it	Istruzioni per l'uso	51
lt	Naudojimo instrukcija	55
lv	Lietošanas instrukcija	59
nl	Gebruiksaanwijzing	63
no	Bruksanvisning	67
pl	Instrukcja stosowania	71
pt	Manual de instruções	75
ro	Instrucţiuni de utilizare	79
ru	Инструкция по применению	83
sk	Návod na použitie	87
sl	Navodila za uporabo	91
sv	Bruksanvisning	95
tr	Kullanım kılavuzu	99

Ordering information

Freeprint® ortho 385 clear		Freeprint® splint 2.0 385 clear	
1000 g	03989	500 g	02080
5 kg Eco Bag	04323	1000 g	02076

Freeprint® denture impact pink-transparent		Freeprint® denture impact pink	
1000 g	04436	1000 g	04437

Freeprint® ibt 385 clear		Freeprint® tray 2.0 385 green	
1000 g	04249	1000 g	02505

**Freeprint® temp 385
500 g**

A1	04058
A2	04059
A3	04060

**Freeprint® temp 385
1000 g**

A1	04062
A2	04063
A3	04064

**Freeprint® splintmaster flex 385
clear**

1000 g	04432
---------------	-------

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Resin for dental 3D printing

Indication

Orthodontic bases, surgical guides and X-ray templates for implant dentistry, occlusal orthos, fixation and transfer keys

Patient Target Group

Persons being treated in the context of a dental procedure.

Intended Users

Dentist, dental technician

Suitable for the following DLP printers/cleaning/post-exposure

see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Bottle: Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Eco Bag: Before use, the material should be shaken intensively and homogenized by tumbling or using a roller mixer with a suitable attachment.
- Minimum material thickness for the design is 1.5 mm
- Polish surface mechanically.
- Processing temperature 23 °C ± 2 °C / 73 °F ± 4 °F.

Freeprint® ortho 385 – Sterilization of Drilling Templates

- Locally applicable regulations and hygiene guidelines are to be observed.
- Important: after main cleaning and prior to post-curing, the components to be sterilized are to be dried, e.g. to be heated in an oven to approx. 40 °C / 104 °F for 30 minutes to remove remaining solvent from the cleaning process
- Freeprint® ortho 385 is validated for the „W&H Lisa 522, Program UNIVERSAL 121“ sterilization process (technical data on the right). Steam sterilization may only be performed with equipment that complies with the standards EN 13060 and EN 285.
- Responsibility for sterility lies with the user. Safety Information

Freeprint® ortho (405 nm)

- Do not use heat-based methods for disinfection or sterilization. This could possibly deform the workpiece.

Technical Data for the sterilisation Process

- A fractionated pre-vacuum with a plateau phase at 121 °C / 250 °F / 15 minutes is to be selected, whereby this refers to the requirements for a B-cycle small steriliser in terms of its venting performance.

Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- MMA-free.
- Refer to the relevant safety data sheet for hazard and safety information.

Cleaning

By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards. Cleaning tablets can be used following the manufacturer's instructions.

Note

- Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- Do not use any acids or solvents.

Note

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Please follow the instructions on the safety data sheet!

For users and/or patients

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- Store Freeprint® ortho in a dry place (at 15 °C - 28 °C / 59 °F - 82 °F) and away from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of Freeprint® ortho may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Freeprint® ortho only insert intraorally in completely polymerised state.

Adverse Effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage

Application

at 23 °C ± 2 °C / 73 °F ± 4 °F

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

① Construction process

Generation of a Print Job complying with machine and material parameters

② Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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detax GmbH