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Light-curing resin based on (meth)acrylate,
for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

Suitable for the following DLP printers/cleaning/post-exposure
see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description).
- Bottle: Before use, the material should be shaken intensively and homogenized with a bottle roller.
- Ecobag: Before use, the material should be shaken intensively and homogenized by kneading.
- Processing temperature $23\text{ °C} \pm 2\text{ °C}$ / $73\text{ °F} \pm 4\text{ °F}$.

Safety Information

- Only to be used by trained specialists for the specified area of application.
- Product contains (meth)acrylates and phosphine oxides. The ingredients of dx model basic may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly closed, carefully close immediately after each use.
- Product may cause allergic reactions.
- Read and understand the safety data sheet!

Storage

- dx model basic is to be stored dry (at $15\text{ °C} - 28\text{ °C}$ / $59\text{ °F} - 82\text{ °F}$) and protected from light. Minimal influence of light can already induce polymerisation.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Area of application

3D Print working and presentation models for the entire spectrum of models in dentistry

Processing

at $23\text{ °C} \pm 2\text{ °C}$ / $73\text{ °F} \pm 4\text{ °F}$

Storage

Ordering information

dx model basic	caramel
1000 g	04734
5 kg Ecobag	04735
dx model basic	light tan
1000 g	04737
5 kg Ecobag	04738
Freeprint® gingiva 385	gingiva
500 g	02820
1000 g	02843

dx model basic

Manufacturing process

① Preparation

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

② Construction process

Generation of a Print Job complying with machine and material parameters

③ Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended.

If possible, post-processing should commence immediately following the construction process.

④ Cleaning

see "Annex 1, Cleaning Equipment"

⑤ Post-exposure

see "Annex 1, Light curing Equipment"