3D Tec Material

rapid cast 2.0

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rapid cast 2.0

UV-curing resin based on (meth) acrylate for the generative fabrication of casts for soft silicone earmolds, for DLP-systems with UV / LED 385 nm

Suitable for the following DLP printers/cleaning/post-exposure see "Annex 1" (enclosed separately)

Processing

- Observe process description.
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.

Safety Information

- Only for the specified intended use by trained specialists.
- Contains (meth)acrylates. Some ingredients of rapid cast 2.0 may cause allergic reactions in predisposed persons. In such cases refrain from using the product.
- Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, safety goggles) when handling the liquid material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

Important notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Product may cause allergic reactions.
- Not suitable for fabricating ear moulds.
- Read and understand the safety data sheet!

Storage

- Rapid cast 2.0 is to be stored dry (at 15 °C - 28 °C / 59 °F - 82 °F) and protected from light. Minimal influence of light can already induce polymerisation.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Disposal

- Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Area of application

Generative fabrication of casts


Ordering information

rapid cast 2.0	1000 g
blau-transparent	04740
detax softwear® 2.0	
Cartridge pack 40 Shore A	8 × 50 ml
clear-transparent	04230
rose-transparent	04231
Cartridge pack 40 Shore A	2 × 50 ml
red opaque	04232
blue opaque	04233
green opaque	04235
yellow-fluorescent	04234
pink-fluorescent	04236
night glow	04229
Cartridge pack 60 Shore A	8 × 50 ml
klar-transparent	03906
rosé-transparent	03911
Cartridge pack 60 Shore A	2 × 50 ml
red opaque	04242
blue opaque	04238
green opaque	04240
yellow-fluorescent	04239
pink-fluorescent	04241
night glow	04237
Mixing tips	25 pcs.
rose	02770
Mixing gun	1 pc.
	04192
softwear coat	500 ml
	03509

Processing

at 23 °C ± 2 °C / 73 °F ± 4 °F

Storage

3D Tec Material

rapid cast 2.0

Manufacturing process for casts Step by step

① Design process

Data preparation according to instructions of the CAD manufacturer.
Wall thickness cast shell: min 0.8 mm
Support: Use a filling cone or optionally a second stand as a support cone.

② Construction process

Generation of a Print Job complying with machine and material parameters
(see „Annex 1, Digital light Processing (DLP) printer, operation software and parameter).

③ Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

④ Cleaning

see „Annex 1, Cleaning Equipment“

⑤ Drying

Leave the components at room temperature for 15 minutes before post curing to remove the solvent residue from cleaning process.

⑥ Post-exposure

see „Annex 1, Light curing Equipment“

⑦ Filling of the cast

Then the cast can be filled with silicone, i.e. detax softwear® 2.0.
Working pressure pneumatic mixing gun: 2 bar.

⑧ Setting process

Put the silicone-filled cast moulds for at least 30 minutes place in a pressure vessel filled with water at 35 °C - 40 °C / 95 °F - 104 °F and 2 - 4 bar. The filled cast moulds must be completely covered with water.

⑨ Demoulding process

After approx. 30 minutes in the warm (35 °C - 40 °C / 95 °F - 104 °F) water, the silicon earmoulds are ready to be removed from the cast moulds. No tools are required.

detax GmbH

Carl-Zeiss-Str. 4 • 76275 Ettlingen • Germany
+49 7243 510 0 • post@detax.com • detax.com