

Workflow Guide · dx rapid cast

3D printing of cast shells that can be removed without any tools
– the dx rapid cast material.



This guide will lead you through the process of designing, manufacturing, preparing, filling and removing the cast shell.

Workflow

1. Design



Time

Approx. 10 min



Design Parameter

**Compatible with all types
of design software (no
additional settings required)**



Wall Thickness Cast

min. 0.8 mm



Filling Cone

**Use a filling cone or second stand
as a support cone.**

Make sure that you place the filling cone on top of the cast shell.

2. Create Print Job



Time

Approx. 5 min



Note

Use the validated parameter files for Rapid Cast.

Please use thicker support bars as by hard cast shells.

3. Printing Process



Time

Approx. 45 min



Note

Use a Detax-validated 3D printer.

4. Post Process

4.1. Cleaning



Time

- **Ultrasonic bath: 2 x 2 min in iso-propyl alcohol (purity $\geq 98\%$)**
Please use two separate containers.
- **Software controlled device:**
Use the following settings:
Detax Rapid cast



Note



- **Between the two cycles and after the washing process, careful blow out with compressed air with max. 2 bar.**
- **Between the two cycles and after the washing process, careful blow out with compressed air with max. 2 bar.**



Attention

Handle the cast molds always with care.

Don't remove the cast shell from the supports, wash it in one part.

Please make sure that the cast shell is complete floated with IPA during the washing process (no air bubbles inside).

4.2. Post Curing



Time

Depending on device

Otoflash G171 N²

2 x 1000 flashes with nitrogen

Rapidshape RScore

Medicalprint cast

Dreve PCU LED

2:30 min with Nitrogen
Light intensity 33%

If necessary, carefully remove the cast shells from the support structure after post-curing.
Don't pull them off, use a cutter.

5. Filling the Cast Molds



Note

No cast separator is required.



Pressure

Working pressure of the pneumatic mixing gun is 2 bar.

6. Setting of the Silicone and Removal Process



Time

- **Softwear 2.0**
- **Approx. 15 min**



Water Temperature

35-40°C



Pressure

2-3 bar

Place the filled cast shells in a pressure pot filled with water to cure. The cast shells must be surrounded by water.

Please check in the IFU of the silicone.

7. Demolding Process

The demolding process starts when the silicone starts to cure.



Time

Few seconds



Tools

Your hands

No tools are necessary to remove the cast shells.

After the specified curing time of the silicone, you can remove the cast shells from the pressure pot and start to remove the cast shell.
