









Date: 07/2025

dx validation printer matrix (385 nm)






















	denture/C&B						splint/surgical guide				model				tray/ibt/cast		
																	
																	
ivoclar																	
PrograPrint PR5				●	●	●	●	●	●	●				●	●		●
MICROLAY <sup>™</sup> 3D-PRINTERS																	
Versus	●			●	●	●	◐	◐	●	●	●			●	●		●
Glidewell																	
Fastprint.io				●													
MICRAFT																	
Prime/Hyper Series	●				●				●	●	●			●	●		●
Ultra Series	●				●				●	●				●	●		
Alpha	●				●		●	●	●	●	●		●	●	●		●
HEYGEARS																	
A2D	●	●	◐	◐	●		◐	◐									
A2D HD		●	◐	●			◐	◐									
A3D		◐	◐	◐			◐	◐									
Chairside		◐	◐	◐	◐	◐	◐	◐									
W2P																	
SolFlex		●			●	●			●	●	●			●	●		●
Qualification	● validated   ◐ in process																

dx validation printer matrix (405 nm)

	denture/C&B						splint/surgical guide				model				tray/ibt/cast		
formlabs																	
Form 4B		●	●	●	●	●	●	●									
RAYSHAPE																	
Edge E2		●	●	●	●		●	●									
MICROLAY																	
Eve Pro											●		●	●			●
Aidite																	
CPD-100	●			●	●		●	●									
phrozen																	
Sonic 4K		●			●					● 405 nm	●	●	●	●	●		●
SHINING 3D DENTAL																	
Accu-Fab L4D		●	●	●	●					● 405 nm	●	●	●	●			●
AccuFab CEL		●	●	●	●	●				● 405 nm	●	●	●	●			●
Accu-Fab D1s		●									●	●	●	●			●
Qualification	● validated   ● in process																



dx validation curing matrix

	denture/C&B						splint/surgical guide				model				tray/ibt/cast		
																	
																	
<b>NK-Optik</b>																	
Otoflash G171 N2	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
NK Flash 250/500	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
NK Flash 150	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐
<b>ASIGA</b>																	
Asiga Cure		◐	◐	◐	◐	◐	◐	◐		◐	◐	◐	◐				◐
<b>rapidshape</b>																	
RS Cure	●	●	●	●	●	●	●	●		●	●	●	●	●	●	●	●
RS Cure XL	●	●	●	●	●	●	●	●		●	●	●	●	●	●	●	●
<b>straumann</b>																	
P Cure	●	●	●	●	●	●	●	●		●	●	●	●	●	●	●	●
<b>dentalfarm</b>																	
Photopol	●				●	●			●		●		●	●	●	●	●
<b>Glidewell</b>																	
ILCD				●													
<b>MTC</b>																	
BB-Cure											●		●	●			●
Qualification	● validated ◐ in process																

dx validation curing matrix

denture/C&B						splint/surgical guide				model				tray/ibt/cast		
phrozen																
Phrozen Cure										●						
ivoclar																
PrograPrint Cure						●								●		
formlabs																
Form Cure (2nd Generation)						●										
Form Cure										●						
HEYGEARS																
PCU Pro						●										
PCU 3.0						●										
Dreve																
PCU LED N2																
SHINING 3D DENTAL																
FabCure										●						
PRUSA RESEARCH																
MedicalCW One										●						
Qualification																