

dx validation printer matrix (385 nm)

	denture/C&B						splint/surgical guide				model					tray/ibt/cast		
ASIGA																		
Max	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Max 2	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Ultra	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Pico2	●				●	●			●	●	●		●	●	●	●	●	●
PRO2	●				●	●			●	●	●		●	●	●	●	●	●
PRO 4K	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
rapidshape																		
ONE		●	●	●	●		●	●									●	
Pro20	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
Pro30	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
D10/D20 Series	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
D30/D40 Series	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
D50+ Series	●	●	●	◐	◐	◐	●	◐	◐	◐				◐	◐	◐	◐	◐
D70/D90 Series					●	●				●				●				
straumann																		
P Series	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
Qualification	● validated ◐ in process																	



dx validation printer matrix (385 nm)

	denture/C&B						splint/surgical guide				model					tray/ibt/cast		
																		
																		
ivoclar																		
PrograPrint PR5				●	●	●	●	●	●	●					●		●	
MICROLAY [®] 3D PRINTERS																		
Versus	●	◐	◐	●	●	●	◐	◐	●	●	●				●	●	●	
Glidewell																		
Fastprint.io				●													◐	
LuxCreo																		
iLux Pro Dental		◐	◐	◐			◐											
MICRAFT																		
Prime/Hyper Series	●				●				●	●	●			●	●	●	●	●
Ultra Series	●				●				●	●						●	●	
Alpha	●				●		●	●	●	●	●	●		●	●	●	●	●
HEYGEARS																		
A2D	●	●	●	◐	●		●	●										
A2D HD		●	●	●			●	◐										
A3D		◐	◐	◐			●	◐										
Chairside		◐	◐	◐	◐	◐	◐	◐										
W2P																		
SolFlex	●	●		●	●	●			●	●	●			●	●	●	●	●
Qualification	● validated ◐ in process																	

dx validation printer matrix (405 nm)

	denture/C&B						splint/surgical guide				model					tray/ibt/cast		
																		
																		
formlabs																		
Form 4B		●	●	●	●	●	●	●				◐	◐					
RAYSHAPE																		
Edge E2		●	◐	●	●		◐	◐										
Edge Mini		◐		●														
MICROLAY																		
Eve Pro											●			●	●			●
Aidite																		
CPD-100	●	◐		●	●		◐	◐										
phrozen																		
Sonic 4K		●								405 nm ●	●	●		●	●	●		●
SHINING 3D																		
Accu-Fab L4D		●	◐	●	●					405 nm ●	●	●		●	●			●
AccuFab CEL		●	◐	●	◐	◐				405 nm ◐	◐	◐		◐	◐			◐
Accu-Fab D1s		◐									●	◐		●	●			●
3D MTC																		
Meccatronicore							◐			405 nm ◐						◐		◐
Qualification	● validated ◐ in process																	

dx validation curing matrix

	denture/C&B						splint/surgical guide				model					tray/ibt/cast		
																		
																		
NK-Optik																		
Otoflash G171 N2	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
NK Flash 250/500	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
NK Flash 150	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐	◐		◐	◐	◐	◐	◐
ASIGA																		
Asiga Cure		●	●	●	●	●	●	◐		●	◐	◐		◐	◐			◐
rapidshape																		
Cure, RSCure, ProCure	●	●	●	●	●	●	●	●		●	●	●		●	●	●	●	●
RS Cure XL	●	●	●	●	●	●	●	●		●	●	●		●	●	●	●	●
straumann																		
P Cure	●	●	●	●	●	●	●	●		●	●	●		●	●	●	●	●
dentalfarm																		
Photopol	●				●	●			●		●			●	●	●	●	●
Glidewell																		
ILCD				●													◐	
LuxCreo																		
ILCP		◐	◐	◐			◐											
Qualification	● validated ◐ in process																	

dx validation curing matrix

	denture/C&B						splint/surgical guide				model					tray/ibt/cast		
																		
																		
																		
BB-Cure											●			●	●			●
BB-Cure N							◐			405 nm ◐						◐		◐
																		
Phrozen Cure											●			●	●			●
ivoclar																		
PrograPrint Cure				●	●	●	●	●	●	●					●	●		●
formlabs 																		
Form Cure (2nd Generation)		◐	◐	◐	◐	◐	◐	◐										
Form Cure											●				●	●		●
																		
PCU Pro		●	●	◐			◐											
PCU 3.0	●	●	◐	●	●		◐											
Dreve																		
PCU LED N2		●	●	●			●	◐										
SHINING 3D DENTAL																		
FabCure											●			●	●			●
Qualification	● validated ◐ in process																	

tray/ibt/cast



● validated ◐ in process

405 nm

Date: 01/2026